

HELIUM/HYDROGEN SEPARATION IN HIGH CONCENTRATION LEVEL

Technique : Micro-GC

Column : Tamis 10m

Carrier gas : Argon


Injector : BF

Column temperature : 45°C

Backflush time : 10 s

Column pressure : 28 psi

μTCD sensitivity : Standard



The Micro-GC is a powerful solution that provides fast, accurate and easy-to-use solution to analyze your gas sample.

A wide range of capillary columns and MEMS injectors are available to optimize channel combinations inside the Micro-GC to obtain the desired separation. The μTCD is the high-performance universal detector used in such systems. A parallel analytical configuration of the channels or analytical modules will allow the sampling of a small amount of gas sample and simultaneous injection in all channels to develop up to 4 chromatograms and one cumulative report.

SRA Instruments is the ideal Partner with specific competences to develop with you turn key solutions based on optimized Micro-GC, accessories and dedicated software utility or report calculations for specific plant context.

#	RT (sec)	Concentration (%)
He	44.7	73.14%
H ₂	47.06	27.78%
O ₂	63.15	0.02%
N ₂	101.62	0.06%



More information : www.sra-instruments.com