


NEW LAB
SOLUTIONS



LARGE VOLUME HEADSPACE THERMAL BATH

PRODUCT FEATURES

- ⬡ Six position tray from 440 mm to 1000 ml jar
- ⬡ Thermal bath for a uniform heat distribution
- ⬡ Temperature set by Chronos software within the sample preparation method
- ⬡ Special jar design for headspace sampling and to have a stable positioning on heated tray
- ⬡ Up to a temperature of + 100 ° C (for higher temperatures, solutions with stainless steel bath can be evaluated)
- ⬡ Customizable bath dimension.

SRA 

INSTRUMENTS
ANALYTICAL SOLUTIONS

www.srainstruments.it

 **CHROMLINE**

SRA Instruments Spa - IT
Ph.: +39 02.92143.58
info@srainstruments.com

SRA Instruments Sas - FR
Tel: +33 4.78442947
info@sra-instruments.com