


NEW LAB  
SOLUTIONS



## LARGE VOLUME HEADSPACE THERMAL BATH

### PRODUCT FEATURES

- ⬡ Six position tray from 440 mm to 1000 ml jar
- ⬡ Thermal bath for a uniform heat distribution
- ⬡ Temperature set by Chronos software within the sample preparation method
- ⬡ Special jar design for headspace sampling and to have a stable positioning on heated tray
- ⬡ Up to a temperature of + 100 ° C (for higher temperatures, solutions with stainless steel bath can be evaluated)
- ⬡ Customizable bath dimension.

**SRA**   
**INSTRUMENTS**  
ANALYTICAL SOLUTIONS

[www.srainstruments.it](http://www.srainstruments.it)

 **CHROMLINE**

SRA Instruments Spa - IT  
Ph.: +39 02.92143.58  
[info@srainstruments.com](mailto:info@srainstruments.com)

SRA Instruments Sas - FR  
Tel: +33 4.78442947  
[info@sra-instruments.com](mailto:info@sra-instruments.com)


NEW LAB  
SOLUTIONS



## SDLE - SPME DUAL LAYER EXTRACTION

### PRODUCT FEATURES

- Easy installation above the standard heated agitator
- Temperature set by Chronos software within the sample preparation method
- Special vial cap to better heat transfer surface
- Up to - 20 °C (this specific are defined by the thermo exchanger installed)



**SRA**   
**INSTRUMENTS**  
ANALYTICAL SOLUTIONS

[www.srainstruments.it](http://www.srainstruments.it)

 **CHROMLINE**

SRA Instruments Spa - IT  
Ph.: +39 02.92143.58  
[info@srainstruments.com](mailto:info@srainstruments.com)

SRA Instruments Sas - FR  
Tel: +33 4.78442947  
[info@sra-instruments.com](mailto:info@sra-instruments.com)